

**Press Release: Solvotrode easyClean – cleaning at the touch of a button!**

Date: 21st September 2010

**The proven Metrohm Solvotrode is now available with the new *easyClean* diaphragm: just press on the electrode head and the electrolyte flows out. The diaphragm does not need to be touched any more for further cleaning. The user benefits.**

- easy contact-free cleaning of the diaphragm
- greater accuracy and reproducibility of the electrolyte flow (glass on glass, the spring in the electrode head returns to its starting position)
- reduced immersion depth of the sensor



Visit [www.metrohm.ie](http://www.metrohm.ie) for details or contact us on 01928 579 600

For more information on this press release please contact the Marketing Department at Metrohm Ireland Ltd.

**Helen Davey - Marketing Manager**  
**Metrohm Ireland Ltd**  
**Tel. 01928 579600**  
**Email [hdavey@metrohm.co.uk](mailto:hdavey@metrohm.co.uk)**

**Note to Editors:**

Metrohm Ireland Ltd is the Irish distributor of Metrohm AG, the leading manufacturer of Ion Analysis Instruments.

**Metrohm** is synonymous with potentiometric titration, Karl Fischer titration, ion chromatography and voltammetry solving your analytical chemistry problems with quality instrumentation and support services in the technologies of POTENTIOMETRIC TITRATION, KARL FISCHER TITRATION (moisture analysis), ION CHROMATOGRAPHY, VOLTAMMETRY (trace analysis), On-line & AT line analysis and VISCOMETRY.

Set up in 1996, Metrohm Ireland serves the key markets of the Ireland including the chemical, petrochemical and pharmaceutical industries, as well as the automotive and food sectors.